

HOW TO USE FENTANYL TEST STRIPS

Fentanyl is the leading cause of overdose in New England and can be found in **ANY** drug.

Test each drug. Every time. Carry Narcan. Never us alone

1 CRUSH AND MIX WHOLE PILL OR ROCK (BEST PRACTICE)

If this is not preferred, bag residue or a piece of a drug can be tested.



2 COMBINE PINCH OF SUBSTANCE WITH 5ML WATER (ONE TUBE) IN COOKER



The amount of drug pictured on the pencil is enough. For residue, place 5mL water (one tube) in baggie.



3 MIX WATER WITH TEST STRIP

Hold the blue end of the test strip and insert the other end into the liquid, no higher than the blue line. Stir for 15 seconds.




4 WAIT UP TO 5 MINUTES FOR RESULTS

Set the strip down on a flat surface and wait.



- **One red line on top is a POSITIVE result** and means the substance likely contains fentanyl.
- **Two red lines is a NEGATIVE result** and means the substance likely does not contain fentanyl.

No red lines means the test is invalid. Be careful: fentanyl test strip results are not 100% accurate.



TESTING DRUGS CONTAINING METH, MDMA, OR ADDERALL

Testing these drugs requires using more water to avoid false positive results.

Follow the steps above but then add a half a cup of water.

